## en Enteral syringe

**Intended use:** KD-JECT®E enteral syringes are intended for the enteral administration of fluids, liquid diets or drugs and for enteral aspiration or irrigation procedures for single use on a single patient only.

### Indications:

- Enteral administration of liquid diets or drugs via feeding tubes/extensions tubes equipped with ENFit® compatible male connectors.
- Enteral aspiration and irrigation procedures via tubes equipped with ENFit® compatible male connectors.

**Intended users:** Qualified healthcare professionals.

**Description of the product:** KD-JECT®E enteral syringes are provided in sterile condition for transient use equipped with a female ENFit® connector; The scale allows exact measuring of the volumes to be administered or aspirated (accuracy of the scale in accordance with ISO 7886-1); KD-JECT®E enteral syringes consist of barrel, purple plunger, plunger stopper and optional protective cap.

**Materials used:** Polypropylene (PP), Latex-free isoprene rubber, silicone oil (lubricant), purple pigment.

### Contraindications:

- Not to be used for parenteral applications.
- Not to be used for patients with known hypersensitivity to any of the materials used.

### Instructions for use:

1. Follow the institutional protocol for using enteral syringes.
2. Select appropriate size of enteral syringe.
3. Connect the syringe to an ENFit® compatible male connector during filling, administration and any other plunger operation to reduce the risk of inaccurate dosage.
4. Fill the syringe by using a bottle adapter or straw – it is not recommended to fill the syringe by diving the syringe tip into the fluid as the tip/the cavity surrounding the tip in the low-dose syringes is to be kept free from any excess fluid.
5. Fully depress the plunger prior to filling and then draw the fluid into the syringe by using the patient-side edge of the black gasket for measuring the volume. Do not fill beyond the maximum scale graduation.
6. Tap the syringe while holding it with the tip pointing up and slightly draw back and push the plunger in order to expel air bubbles. Repeat until the air is fully ejected.
7. If there is excess fluid in the tip/the cavity surrounding the tip in the low-dose syringes, tap the syringe to completely remove any excess fluid.
8. Use the syringe immediately after preparation.

### ⚠ Cautions:

- Use of this product is restricted to qualified healthcare professionals.
- Read instructions before use.
- Visually inspect and check the expiry date as well as the product and package integrity prior to use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- Do not use if the syringe is damaged/deformed/blocked.
- Do not use if the individual package is damaged.
- Do not use after the expiry date.
- Use the product immediately after opening the individual package.
- For single use on a single patient only. Re-use of single use devices creates a potential risk for patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Discard after single use. Do not recycle, reprocess or re-sterilize.

**Storage conditions:** Store in a clean and dry place. Do not expose to heat or direct sunlight. Do not store at extreme temperatures.

**Disposal:** Discard the product after single use in accordance with pertinent national and institutional guidelines.

**Notice:** Any serious incident that has occurred in relation to the device must be reported to KD Medical GmbH Hospital Products and to the competent authority of the EU Member State (or state with the same rules (Regulation (EU) 2017/745 on medical devices)) in which the user and/or patient is established.

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