en Oral syringe

Intended purpose: KD-JECT®O oral syringes are intended for the oral administration of fluids, liquid diets or drugs for single use on a single patient only.

Indications: KD-JECT®O oral syringes are provided in sterile condition for the oral administration of fluids, liquid diets or drugs for transient use in clinical settings in all patient groups.

Intended users: Qualified healthcare professionals.

Description of the product: KD-JECT®O oral syringes are amber colored for light protection and consist of barrel, plunger, plunger stopper and protective cap with different volumes (1, 3, 5, 10, 20ml). The protective cap might be used to cap the syringe to prevent leakage during transportation. KD-JECT®O oral syringes are equipped with a connector incompatible to Luer connectors, thereby preventing accidental parenteral administration. The scale allows exact measuring of the volumes to be administered (accuracy of the scale in accordance with ISO 7886-1).

Nominal capacity of syringe V ml	Accuracy (Tolerance on any graduated capacity)		Precision (Coefficient of variation CV)	
	Less than half nominal capacity	Equal to or greater than half nominal capacity	Less than half nominal capacity	Equal to or greater than half nominal capacity
V<2	±(1,5% of V + 2% of expelled volume)	±5% of expelled volume	CV < 4%	CV < 2%
2≤V<5	±(1,5% of V + 2% of expelled volume)	±5% of expelled volume		
5≤V<10	±(1,5% of V + 1% of expelled volume)	±4% of expelled volume		
10≤V<20	±(1,5% of V + 1% of expelled volume)	±4% of expelled volume		
20≤V<30	±(1,5% of V + 1% of expelled volume)	±4% of expelled volume		

Materials used:

Polypropylene (PP), Latex-free silicone rubber; silicone oil, amber pigment.

Contraindications:

- Not to be used for parenteral applications.
- Not to be used in unconscious patients and patients unable to swallow.
- Not to be used for patients with known hypersensitivity to any of the materials used.

Instructions for use:

1. Follow the institutional protocol for using oral syringes.
2. Select appropriate sizes of oral syringe for the intended procedure.
3. Fully depress the plunger prior to filling and then draw the fluid into the syringe by using the patient-side edge of the black plunger stopper for measuring the volume. Do not fill beyond the maximum scale graduation.
4. Tap the syringe while holding it with the tip pointing up and slightly draw back and push the plunger in order to expel air bubbles. Repeat until the air is fully ejected.
5. Use the syringe immediately after preparation.

⚠ Cautions:

- Use of this product is restricted to qualified healthcare professionals.
- Read instructions before use.
- Visually inspect and check the expiry date as well as the product and package integrity prior to use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- Do not use if the syringe is damaged/deformed/blocked as this may lead to inadequate medicine delivery.
- Do not use if the product or individual package is impaired, damaged, opened or contaminated as this may lead to contamination and/or impairment of functional capability.
- Do not use after the expiry date.
- Use the product immediately after opening the individual package.
- For single use on a single patient only. Re-use of single use devices creates a potential risk for patient or user. It may lead to contamination and/or impairment of functional capability.
- Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Discard after single use. Do not recycle, reprocess or re-sterilize.

Storage conditions: Store in a clean and dry place. Do not expose to heat or direct sunlight. Do not store at extreme temperatures.

Disposal: Discard the product after single use in accordance with pertinent national and institutional guidelines.

Notice: Any serious incident that has occurred in relation to the device must be reported to KD Medical GmbH Hospital Products and to the competent authority of the EU Member State (or state with the same rules (Regulation (EU) 2017/745 on medical devices)) in which the user and/or patient is established.

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