
en Hypodermic needle with needle protection device

Brief description of the product: The safety mechanism of KD-FINE® SAFETY covers the needle after use. In the activated position, the safety shield guards against accidental needlestick injuries during handling of the used product and its disposal.

Materials used: Stainless steel, polypropylene (PP), silicone oil (lubricant), epoxy (adhesive).

Intended purpose: KD-FINE® SAFETY hypodermic needles are medical devices intended for tissue penetration in order to provide a channel for injections or blood, body liquid or cell collection. KD-FINE® SAFETY hypodermic needles are equipped with a female Luer connector and are to be used in conjunction with devices equipped with a male Luer connector in compliance with ISO 80369-7 (such as syringes).


Indications: Intradermal, subcutaneous, intramuscular, intravenous, intraarticular and other soft tissue injections; blood / body liquid / cell collection.

Contraindications: -

Intended user: The use of this product is restricted to qualified healthcare professionals.


Intended patient population: There is no particular intended patient population for the use of Hypodermic Needles. The different sizes and variants available are to be selected by experienced healthcare professionals regarding suitability for the procedure envisaged and the patient's body condition.

Instructions for use:


1
Carefully remove the safety needle from the package and securely attach it to a syringe with a push and clockwise twisting motion.


When used with a Luer slip syringe, make sure the needle-syringe assembly is firm.


2
Move the safety shield away from the needle and towards the syringe barrel. The safety shield will remain in the position you set.

3
Remove the needle cover carefully to avoid damaging the needle. Remove any air before use. Precisely adjust the medication by pushing the plunger until a medication drop emerges at the needle tip. Perform the intended procedure in accordance with currently accepted medical techniques.



4
Immediately after completion of the procedure position the safety shield in preparation for device activation: Using a one-handed technique, push the shield forward with thumb or finger so that the shield is less than 90° from the needle. **Activation can be made in two different ways (A or B):**


A
Center the thumb on the finger pad of the safety shield (5a). Activate the safety shield: Press down with a firm, quick motion until a distinct audible click is heard (6a).

NOTE: The activation force shall be applied vertically in the activation direction of the safety shield.

The device may detach when using Luer slip syringe if improper or wrong activation procedure is used.

B
Position the safety shield in the direction of activation at the lowest possible angle to a flat surface (5b). Activate the safety shield: Press down with a firm, quick motion until a distinct audible click is heard (6b).

The audible click may be missed for small needle sizes. Visually confirm that the needle is fully engaged under the lock (7). *The location of the lock differs according to the needle size. **NOTE:** If the activation failed, dispose of the device immediately into an approved sharps container.

8. Dispose of the used product following federal and local regulations for sharps disposal.

⚠ Precautions, warnings:

- Not validated for intraocular use.
- Read instructions before use. KD-FINE® SAFETY should be used according to these instructions for use.
- KD-FINE® SAFETY is non-toxic, sterile and non-pyrogenic within shelf life if the individual package has not been opened or damaged.
- Visually inspect and carefully check the device and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- Do not use if the product or individual package is impaired, damaged, opened or contaminated.
- The product should be used immediately after opening of the package.
- Do not bend the needle. If the needle is bent or damaged, do not try to straighten or use the product.
- Handle with care to avoid needlestick injuries. Keep hands behind the needle at all times during use and disposal.
- Do not replace the original needle cover after use.

- Activation of the safety mechanism may cause minimal splatter of any fluid remaining on the needle after injection. For maximum safety, activate the safety mechanism by using a one-handed technique away from self and others.
- Do not attempt to deactivate the safety mechanism by forcing the needle out of the safety shield.
- For single use on a single patient only.
- Re-use of single-use devices creates a potential risk for patient or user. It may lead to contamination and/or impairment of functional capability.
- Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Dispose of after single use following federal and local regulations for sharps disposal by using an appropriate container for sharp medical waste.
- Do not clean or re-sterilize.
- KD-FINE® SAFETY must not be re-used.
- Store in a cool and dry place.
- Do not expose to heat or direct sunlight.

Notice: Any serious incident that has occurred in relation to the device must be reported to KD Medical GmbH Hospital Products and to the competent authority of the EU Member State (or state with the same rules (Regulation (EU) 2017/745 on medical devices)) in which the user and/or patient is established.

	Manufacturer
	Catalogue number
	Batch code
	Date of manufacture
	Use-by date

	Fragile, handle with care
	Keep dry
	Keep away from sunlight
	Do not use if package is damaged
	Do not re-use

	Do not resterilize
	Consult instructions for use
	Caution
	Latex-free
	Non-pyrogenic

	Sterilized using ethylene oxide
	Medical device
	Single sterile barrier system
	Unique device identifier
	CE marking

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